



Easy to use


Penetrating trauma




Bleeding starts




Step 1 - Insert applicator




Step 2 - Push in plunger to dispense Celox




Step 3 - Withdraw applicator and cover




Step 4 - Compress



Stable Wound



Step 5 - Wrap for evacuation



Medtrade Products Ltd.  
Crewe Business Park, Crewe, CW1 6GL UK  
T: +44 (0)1270 500019 F: +44 (0)1270 500045

Contact us at [www.celoxmedical.com](http://www.celoxmedical.com)

